



REVOLUTION WIRELESS PROGRAMS AND PROTOCOL PARAMETERS

Using electrotherapy to facilitate movement, strengthen muscle, and improve balance and coordination.

- Wireless allows for the utilization of electrotherapy during functional activity without the restriction of lead wires
- Utilizes VMS and VMSFR
- Up to 4 channels that can run concurrently
- Electrodes have a unique snap
- Modules snap onto the electrodes
- Modules allow for adjustment of lead wire length
- Wraps are utilized to hold the modules in place during activity

Pain:

Acute Pain (Gate control)

Chronic Pain (Opiate Release)

Trigger Point

Edema:

Muscle Pump Continuous

Muscle Pump Cycled

Relax Muscle Spasm/ROM:

Decrease Muscle Tone

Increase ROM

Muscle Strengthening:

Atrophy

Slow Twitch

Fast Twitch

Muscle Re-Education:

VMSFR Static 4 Channel

VMSFR Static 2 Channel

VMSFR Dynamic 4 Channel

VMSFR Dynamic 2 Channel

Pain Management

- Acute Pain
Single Ch; PD=200 μ s; CC; Freq=100 pps; Continuous; 30 mins
- Chronic Pain
Single Ch; PD=300 μ s; CC; Freq=6 pps; Continuous; 30 mins
- Trigger Point
Single Ch; PD=400 μ s; CC; Freq=6 pps; Continuous; 4 mins

DJO, LLC
1430 Decision Street
Vista, CA 92081-8533

760.727.1280
800.321.9549 toll free
www.djortho.com



NEVER STOP GETTING BETTER

Decrease Edema

- Muscle Pump Cycle
Single Ch; PD=200 μ s; CC; Freq=35 pps; Cycle time 5/5; additional Ramp
up/down=2 sec; 20 mins
- Muscle Pump Continuous
Single Ch; PD=200 μ s; CC; Freq=5 pps; Cycle time=Continuous; 20 mins

Relax Muscle Spasm/ROM

- Decrease Muscle Tone
Single Ch; PD=300 μ s; CC; Cycle Time=5:5; Freq=100 pps; additional Ramp
up/down=2 sec; 30 mins
- Increase ROM
Single Ch; PD=200 μ s; CC; Cycle Time=4:12; Freq=35 pps; additional Ramp
up/down=2 sec; 20 mins

Muscle Strengthening

- Muscle Atrophy
Single Ch; PD=300 μ s; CC; Cycle Time=10:50; Freq=50 pps; additional Ramp
up/down=2 sec; 20 mins
- Slow Twitch function
Single Ch; PD=300 μ s; CC; Cycle Time=4:12; Freq=35 pps; additional Ramp
up/down=2 sec; 20 mins
- Fast Twitch function
Single Ch; PD=200 μ s; CC; Cycle Time=5:5; Freq=50 pps; additional Ramp
up/down=2 sec; 20 mins

Muscle Re-Education

- VMS-FR Static 2 CH.
Ch 1+2; PD=300 μ s; CC; Burst Duration=4500 msec; On/Off time=1:3; Freq=35
pps; 30 mins
- VMS-FR Static 4 CH.
Ch 1+2; PD=300 μ s; CC; Burst Duration=4500 msec; On/Off time=1:3; Freq=35
pps; 30 mins
- VMS-FR Dynamic 2 CH.
Ch 1+2; PD=200 μ s; CC; Burst Duration=2000 msec; On/Off time=1:5; Freq=80
pps; 20 mins
- VMS-FR Dynamic 4 CH.
Ch 1+2; PD=200 μ s; CC; Burst Duration=2000 msec; On/Off time=1:5; Freq=80
pps; 20 mins

DJO, LLC
1430 Decision Street
Vista, CA 92081-8533

760.727.1280
800.321.9549 toll free
www.djortho.com

DECREASE PAIN

IFC Gate 40

Beat Freq. 80-150 bps
Scan 40%
Carrier Freq. 4000
CV

IFC Gate 100 (ACP IFC Sensory)

Beat Freq. 80-150 bps
Scan 100%
Carrier Freq. 4000
CV

IFC Opiate 40

Beat Freq. 1-10 bps
Scan 40%
Carrier Freq. 2500
CV

IFC Opiate 100 (ACP IFC Motor)

Beat Freq. 1-10 bps
Scan 100%
Carrier Freq. 2500
CV

IFC Nerve Block 100 (ACP IFC Nerve Block)

Sweep OFF
Beat Freq. 100 bps
Scan 100%
Carrier Freq. 5000 Genisys 10,000 NEO
CV

IFC Sensory/Motor (ACP Sensory/Motor) Use Sequencing Program

Sensory

Sweep OFF
Beat Freq. 100 bps
Scan 100%
Carrier Freq. 4,000
CV
Treatment Time 15:00

Motor

Sweep OFF
Beat Freq. 2 bps

Vector Scan 100%
Carrier Freq. 5,000 Hz
CV
Treatment Time 15:00

IFC Motor/Sensory (ACP Motor Sensory) Use Sequencing Program

Motor
Sweep OFF
Beat Freq. 15 bps
Vector Scan 100%
Carrier Frequency 5,000 Hz
CV
Treatment Time 8:00
Motor
Sweep OFF
Beat Freq. 2 bps
Vector Scan:Off
Carrier Freq. 5,000 Hz
CV
Treatment Time 7:00

Sensory
Sweep OFF
Beat Freq. 100 bps
Vector Scan 100
Carrier Freq. 4,000 Hz
CV
Treatment Time 15:00

Premod Gate

Freq. 80-150 bps
Carrier Freq. 4,000 Hz (NEO only)
Cycle Time Continuous
CV

Premod Opiate

Freq. 1-10 bps
Carrier Freq. 2500 Hz (NEO only)
Cycle Time Continuous
CV

MUSCLE FUNCTION/CONTROL

VMS-FR Tonic (Weakness, Flaccidity, Abnormal Tone) (ACP PENS)
Phase Duration 200 μ sec (increase for flaccid muscle and larger muscle)
Freq. 50 pps
Burst Duration 2000 μ sec (Set to equal the duration of the functional activity)

On/Off Ratio 1:3
CC
Set Channels Simultaneously

VMS-FR Phasic (Weakness, Dynamic Instability, Poor Power) (ACP PENS)
Phase Duration 200 μ sec (increase for flaccid muscle and larger muscle)
Freq. 80 pps
Burst Duration 2000 μ sec (Set to equal the duration of the functional activity)
On/Off Ratio 1:5
CC
Set Channels Simultaneously

MUSCLE STRENGTHENING

VMS-Power Single
Channel Mode Single
Phase Duration 200 μ sec (increase for flaccid muscle and larger muscle)
Cycle Time 10/50 (cycle time will change determined by patient abilities)
CC
Freq. 50 pps
Anti-Fatigue Off
Ramp 2 sec.

VMS-Power Reciprocal
Channel Mode Reciprocal
Phase Duration 200 μ sec (increase for flaccid muscle and larger muscle)
Cycle Time 10/50 (cycle time will change determined by patient abilities)
CC
Freq. 50 pps
Anti-Fatigue Off
Ramp 2 sec.
Set Channels Independently

VMS-Power Co-Contract
Channel Mode Co-Contract
Phase Duration 200 μ sec (increase for flaccid muscle and larger muscle)
Cycle Time 10/50 (cycle time will change determined by patient abilities)
CC
Freq. 50 pps
Anti-Fatigue Off
Ramp 2 sec.
Set Channels Independently

VMS-Endurance Single

Channel Mode Single

Phase Duration 200 μ sec (increase for flaccid muscle and larger muscle)

Cycle Time 10/50 (cycle time will change determined by patient abilities)

CC

Freq. 35 pps

Anti-Fatigue Off

Ramp 2 sec.

VMS-Endurance Reciprocal

Channel Mode Reciprocal

Phase Duration 200 μ sec (increase for flaccid muscle and larger muscle)

Cycle Time 10/50 (cycle time will change determined by patient abilities)

CC

Freq. 35 pps

Anti-Fatigue Off

Ramp 2 sec.

Set Channels Independently

VMS-Endurance Co-Contract

Channel Mode Co-Contract

Phase Duration 200 μ sec (Increase this as the size of the muscle increases)

Cycle Time 10/50 (increase for flaccid muscle and larger muscle)

CC

Freq. 35 pps

Anti-Fatigue Off

Ramp 2 sec.

Set Channels Independently

Anti-Fatigue Option

Anti-Fatigue ON (frequency modulation starting with a higher frequency which is then reduced over time to prevent fatigue of the muscle fiber and allow for longer treatment times to build endurance)

VMS Burst -Power Single (Better than ACP MFAC)

Channel Mode Single

Phase Duration 200 μ sec (increase for flaccid muscle and larger muscle)

Cycle Time 10/50 (cycle time will change determined by patient abilities)

CC

Freq. 50 pps

Anti-Fatigue Off

Ramp 2 sec.

VMS Burst -Power Reciprocal

Channel Mode Reciprocal

Phase Duration 200 μ sec (increase for flaccid muscle and larger muscle)

Cycle Time 10/50 (cycle time will change determined by patient abilities)

CC

Freq. 50 pps

Anti-Fatigue Off

Ramp 2 sec.

Set Channels Independently

VMS Burst -Power Co-Contract

Channel Mode Co-Contract

Phase Duration 200 μ sec (increase for flaccid muscle and larger muscle)

Cycle Time 10/50 (cycle time will change determined by patient abilities)

CC

Freq. 50 pps

Anti-Fatigue Off

Ramp 2 sec.

Set Channels Independently

VMS Burst -Endurance Single (Better than ACP MFAC)

Channel Mode Single

Phase Duration 200 μ sec (Increase this as the size of the muscle increases)

Cycle Time 10/50 (increase for flaccid muscle and larger muscle)

CC

Freq. 35 pps

Anti-Fatigue Off

Ramp 2 sec.

VMS Burst -Endurance Reciprocal

Channel Mode Reciprocal

Phase Duration 200 μ sec (Increase this as the size of the muscle increases)

Cycle Time 10/50 (increase for flaccid muscle and larger muscle)

CC

Freq. 35 pps

Anti-Fatigue Off

Ramp 2 sec.

Set Channels Independently

VMS Burst-Endurance Co-Contract

Channel Mode Co-Contract
Phase Duration 200 μ sec (increase for flaccid muscle and larger muscle)
Cycle Time 10/50 (cycle time will change determined by patient abilities)
CC
Freq. 35 pps
Anti-Fatigue Off
Ramp 2 sec.
Set Channels Independently

Anti-Fatigue Option

Anti-Fatigue ON (frequency modulation starting with a higher frequency which is then reduced over time to prevent fatigue of the muscle fiber and allow for longer treatment times to build endurance)

Russian-Power Single (ACP MFAC)

Channel Mode Single
Duty Cycle 50% (may be reduced to make Russian slightly more comfortable, however it will reduce total power)

CC
Anti-Fatigue OFF
Cycle Time 4/12 (cycle time will change determined by patient abilities)
Burst Freq. 50 bps
Ramp Time .5 sec

Russian-Power Reciprocal

Channel Mode Reciprocal
Duty Cycle 50% (may be reduced to make Russian slightly more comfortable, however it will reduce total power)

CC
Anti-Fatigue OFF
Cycle Time 4/12 (cycle time will change determined by patient abilities)
Burst Freq. 50 bps
Ramp Time .5 sec
Set Channel Amplitude Independently

Russian-Power Co-Contract

Channel Mode Single
Duty Cycle 50% (may be reduced to make Russian slightly more comfortable, however it will reduce total power)

CC
Anti-Fatigue OFF
Cycle Time 4/12 (cycle time will change determined by patient abilities)

Burst Freq. 50 bps
Ramp Time .5 sec
Set Channel Amplitude Independently

Russian- Endurance Single

Channel Mode Reciprocal
Duty Cycle 50% (may be reduced to make Russian slightly more comfortable, however it will
reduce total power)

CC

Anti-Fatigue OFF

Cycle Time 12/50 (cycle time will change determined by patient abilities)

Burst Freq. 35 bps

Ramp Time 2 sec

Russian-Endurance Reciprocal

Channel Mode Reciprocal

Duty Cycle 50% (may be reduced to make Russian slightly more comfortable, however it will
reduce total power)

CC

Anti-Fatigue OFF

Cycle Time 4/12 (cycle time will change determined by patient abilities)

Burst Freq. 35 bps

Ramp Time 2 sec

Set Channel Amplitude Independently

Russian-Endurance Co-Contract

Channel Mode Reciprocal

Duty Cycle 50% (may be reduced to make Russian slightly more comfortable, however it will
reduce total power)

CC

Anti-Fatigue OFF

Cycle Time 4/12 (cycle time will change determined by patient abilities)

Burst Freq. 35 bps

Ramp Time 2 sec

Set Channel Amplitude Independently

Anti-Fatigue Option

**Anti-Fatigue ON (frequency modulation starting with a higher frequency which is then
reduced over time to prevent fatigue of the muscle fiber and allow for longer treatment
times to build endurance)**

EDEMA CONTROL

VMS Edema Muscle Pump Cycled

Channel Mode Single
Phase Duration 200 μ sec
Cycle Time 10/10
CC
Freq. 35 pps
Anti-Fatigue Off
Ramp 2 sec.

VMS Edema Muscle Pump Continuous

Channel Mode Co-Contract (set channels independently)
Phase Duration 200 μ sec
Cycle Time Continuous
CC
Freq. 5 pps
Anti-Fatigue Off
Ramp 2 sec.

High Volt Edema (ACP Edema)

Polarity Negative
Sweep Continuous
Ramp 2 sec.
Freq. 100 pps
Cycle Time Continuous

HIGH VOLT

High Volt Spasm

Polarity Negative
Sweep Continuous
Ramp 2 sec.
Freq. 100 pps
Cycle Time Continuous

High Volt Increase Circulation (ACP Increase Circulation)

Polarity Negative
Sweep Continuous
Ramp 2 sec.
Freq. 100 pps
Cycle Time Continuous